

Manipulation N°4

A redox titration

Redox titration

1-Introduction :

The titration consists in determining the normality of a reducing solution knowing that of the oxidizing solution. It is proposed to study the oxidation of the Fe^{2+} ion by the permanganate ion MnO_4^- in an acid medium.

This assay is called manganimetry. The oxidizing properties of the permanganate ion are at the origin of manganimetry. The oxidizing form MnO_4^- is purple, the reducing form Mn^{2+} is colorless, which makes it possible to determine the equivalent point without using colored indicators.

2- Goal of the Practical Work:

This involves determining the normality of a solution (FeSO_4), using a solution of Potassium Permanganate (KMnO_4) prepared in the laboratory.


3- Definitions:

3-1. Oxidation: An oxidation is a reaction during which a reactant yields (loses) one or more electrons.




In this reaction, the Fe^{2+} ion is the oxidized form and Fe is the reduced form of the oxidation-reduction couple (or redox couple) (Fe^{2+}/Fe).


3-2. Reduction: A reduction is a reaction in which a reactant captures (takes) one or more electrons.





In this reaction, the metallic copper is the reduced form and the Cu^{2+} ion is the oxidized form of the redox couple (Cu^{2+}/Cu).

3-3. The oxidant: is the reactant capable of causing an oxidation, that is to say that it captures one or more electrons from another reactant.





3-4. The reducer: is the reagent capable of causing a reduction, that is to say that it yields one or more electrons to another reagent.



3-5. Oxidation-reduction couple: Oxidation-reduction couple or redox couple is a couple formed by an oxidant and its conjugated reducing agent.



Examples of redox couple:

redox couple	oxidant	+ n e ⁻	↔	reducing
Cu²⁺/Cu	Cu²⁺	+ 2 e⁻	↔	Cu
Fe²⁺/Fe	Fe²⁺	+ 2 e⁻	↔	Fe
H⁺/H₂	2H⁺	+ 2 e⁻	↔	H₂
Ag⁺/Ag	Ag⁺	+ 1 e⁻	↔	Ag

3-6. Potassium permanganate:

Potassium permanganate (KMnO₄) is a particularly powerful oxidant. It comes in the form of purple crystals composed of potassium ions, K⁺, and permanganate ions, [MnO₄]⁻. It is odorless and has a bitter taste.

In the laboratory, potassium permanganate is used to perform titrations. At equivalence, the solution indeed changes color, going from purple to pink. Potassium permanganate is also used in water treatment since it oxidizes the iron and manganese contained in groundwater. In everyday life, it can be used, in dilute solution, to eliminate black traces left by fungi between bathroom tiles.

4- Experimental part:

4-1. Dosage of oxalic acid by potassium permanganate:

Equipment: Burette, Erlenmeyer, test tube, wash bottle, KMnO₄ solution (0.1N), H₂C₂O₄.2H₂O solution, Distilled water, H₂SO₄ (10%).

- 1- Rinse the material. Burette, Erlenmeyer.
- 2- Fill the burette with the Normality KMnO₄ solution (N_A= 0.1N).
- 3- Take a 250 mL Erlenmeyer flask and put 50 mL of distilled water in it.
- 4- Add (V_B=10 ml) solution (H₂C₂O₄. 2H₂O).
- 5- Add about 5ml of H₂SO₄ at (10%), and heat to 60-70°C,

- 6- Perform a rapid assay to estimate the equivalence volume.
- 7- Note the volume V_A of $(KMnO_4)$ poured.
- 8- Carry out 2 dosing tests.

4-2. Determination of iron in ferrous sulphate by potassium permanganate:

Equipment: Burette, Erlenmeyer, test tube, wash bottle, $KMnO_4$ solution (0.1N), $FeSO_4 \cdot 7H_2O$ solution, Distilled water, H_2SO_4 (10%).

- 1- Rinse the material. Burette, Erlenmeyer.
- 2- Fill the burette with the Normality $KMnO_4$ solution ($N_A = 0.1N$).
- 3- Take a 250 mL Erlenmeyer flask and put 50 mL of distilled water in it.
- 4- Add ($V_B = 10$ ml) solution ($FeSO_4 \cdot 7H_2O$).
- 5- Add about 5ml of H_2SO_4 at (10%), and heat to 60-70°C,
- 6- Perform a rapid assay to estimate the equivalence volume.

5- Results and calculations:

5-1. Dosage of oxalic acid by potassium permanganate:

- 1- Purpose of the practical work.
- 2- Write the oxidation-reduction half-reactions, specify the redox couples.
- 3- Write the overall reaction.
- 4- Note the volume (V_A) of $(KMnO_4)$ poured and calculate the Normality (N_B) of $(H_2C_2O_4 \cdot 2H_2O)$.
- 5- Calculate the Concentration (C_B) of the solution $(H_2C_2O_4 \cdot 2H_2O)$.
- 6- What to conclude when the role of sulfuric acid H_2SO_4 ?
- 7- For what reason in heating.
- 8- Determine the limiting reagent of this reaction.
- 9- Can sulfuric acid be replaced by HCl or H_3PO_4 . Explain.

5-2. Determination of iron in ferrous sulphate by potassium permanganate:

- 1- Purpose of the practical work.
- 2- Write the oxidation-reduction half-reactions, specify the redox couples.
- 3- Write the overall reaction.
- 4- Note the volume (V_A) of $(KMnO_4)$ poured and calculate the Normality (N_B) of $(FeSO_4 \cdot 7H_2O)$.
- 5- Calculate the Concentration (C_B) of the solution $(FeSO_4 \cdot 7H_2O)$.
- 6- Determine the limiting reagent of this reaction.
- 7- Can sulfuric acid be replaced by HCl or H_3PO_4 . Explain.