TP N°5: Titration of vinegar

MANIPULATION N°5

Titration of vinegar

Acid/base titration (Application to the titration of commercial vinegar)

1 - Introduction :

Vinegar is a wine (alcohol) made sour by acetic fermentation and used as a condiment or preservative. Vinegar is considered to be an aqueous solution of acetic acid with a density (d) of about 1. The objective of this practical work is to determine the degree of acidity D or (°) of a vinegar, defined as the mass of pure acetic or ethanoic acid contained in 100g of solution. To do this, the acetic (ethanoic) acid contained in a known volume of vinegar will be determined using a strong base solution of known concentration: sodium hydroxide solution (NaOH).

During titration, the concentration of the titrating solution (in this case NaOH) must be known precisely. In most cases, however, the concentrations of the solutions are not strictly accurate because:

- the solution is prepared from a solid product, the purity of which is not guaranteed by the manufacturer.
- the solution is prepared from a product whose sampling (weighing, volume) cannot be precise.
- the chemical composition of the solution changes over time.

2- The degree of acidity of the vinegar (D) or (°):

The degree of a vinegar **D** is the same number as the **mass**, in **grams**, of **pure** ethanoic acid contained in **100g** of **solution**.

The degree of acidity of the commercialised vinegar is indicated on each bottle. 1% acidity corresponds to 1g of ethanoic acid in 100g of vinegar solution. The percentage of acidity therefore corresponds to a percentage by mass.

For a volume of 100g of vinegar: (Ac: acetic acid)

$$m_{Ac} = n_{Ac} \times M_{Ac} = C_a \times V_a \times M_{Ac}$$
 from the formula $n = C \times V$

For vinegar we also have :
$$\rho = m/V$$
 so $V = m/\rho$: $m_{Ac} = \frac{m}{\rho} \times C \times M_{Ac} = \frac{100}{\rho} \times C \times M_{Ac}$ (g)

Degree of acidity =
$$D^{\circ} = m_{Ac} = \frac{100}{\rho} \times C \times M_{Ac}$$

Example: vinegar at 6° contains 6g of pure acetic (ethanoic) acid per 100g of vinegar solution.

3- Objective of the practical work (TP):

- How to control the quality of a product by titration.
- Products tested: different types of commercial vinegar.



TP N°5: Titration of vinegar

4- Experimental part :

• Titration of commercial vinegar (acetic acid) :

Materials: Burette, Erlenmeyer flask, beaker, Graduated cylinder, wash bottle, distilled water, NaOH solution (0.1N), different types of commercial vinegar (CH₃-COOH) and colour indicators.

1 - Preparation of the dilute vinegar solution:

- 1 Using a pipette, take precisely 10 mL of commercial white vinegar.
- 2- Put it in a 100 mL volumetric flask.
- **3** Complete with distilled water to the mark.
- 4- Cap the flask and agitate to homogenise the solution.

2- Titration in the presence of phenophthalein :

- 1 Rinse the equipment. Burette, Erlenmeyer flask, etc.
- 2- Fill the burette with the NaOH solution of concentration ($N_B = 0.1 N$).
- 3- Take (V_A =10 mI) of the (CH_3 -COOH) solution and place it in a 100 mL Erlenmeyer flask.
- 4- Add about 20 mL of distilled water.
- 5- Add 2 drops of the colour indicator: phenolphthalein.
- 6- Make a rapid titration to estimate the volume of equivalence (V_B) .
- 7- Record the volume V_B of (NaOH) added.
- 8- Using a pH paper or pH meter, read the pH of the solution every 2ml until coloured.
- 9- Make 2 titration tests.

 $M_{CH3-COOH} = 60 \text{ g.mol}^{-1}$, $\rho = 1,01 \text{ kg.L}^{-1}$

3-Résultats et calculs :

- 1 Write down the reaction equation:
- 2- Note the volume V_B of (NaOH) added and determine the concentration C_A of (CH₃-COOH).
- **3-** Calculate the initial molar concentration C_0 of the ethanoic acid (commercial vinegar).
- **4-** Calculate the degree of acidity of the vinegar (D) and compare this result with that on the sticker on the bottle of commercial vinegar.
- 5- Draw the graph representing the value of pH= $f(V_B)$ on millimetric paper.
- **6-** Determine the volume V_E and the **pH** at equivalence using the graph.
- 7- Which of the coloured indicators in the table below is best adapted for identifying equivalence. **Justify** your answer.

Indicator	Acidic colour	Transition zone	Basic colour
Helianthine	Red	3,1 - 4,4	Yellow
Methyl red	Red	4,2 - 6,2	Yellow
Bromothymol blue	Yellow	6,0 - 7,6	Blue
Phenolphthalein	Uncoloured	8,2 - 10,0	Pink-violet

8- Why is distilled water **added**? Does the addition of distilled water **change** the volume V_E added at equivalence?



TP N°5: Titration of vinegar

Write down the reaction	equation :		
	•		
Note the volume V _B of (N volume of (NaOH)?:(V concentration C _A ?:	' _B =)	ine the concentration \mathcal{C}_{A} o	f (CH₃-COOH).
Calculate the initial molar	concentration C 0 of the	ethanoic acid (commercial	vinegar).
Calculate the degree of a	cidity of the vinegar (D)		
			/
Draw the graph represent Determine the volume V _E , pH) = (Which of the coloured inc your answer. Indicator Helianthine Methyl red	and the pH at equivalend		tifying equivalence. Jus Basic colour Yellow Yellow Yellow
Bromothymol blue	Yellow	6,0 - 7,6	
Phenolphthalein		0,0 - 7,0	Blue
	Uncoloured	8,2 - 10,0	Blue Pink-violet
Why is distilled water ad	Uncoloured	hanged the volume VE pour	Pink-violet red at equivalence?
Why is distilled water ad	Uncoloured	hanged the volume V _E pour Name:	Pink-violet red at equivalence?





TP N°5: Titration of vinegar

Draw the graph representing the value of $pH=f(V_B)$ on millimetric paper.

													11111111111	
						 111111111	 					1111111111	110000000	
1						 	 		11111111	 		*****		
		:::::::::::::::::::::::::::::::::::::::				 	 			 				
												11111		

	1771													
4 1 4 1						 								
1577				17777		 	 	 *********			1711			
		117711777				 11111	 	 		 				
1-1-1-1					*****		 	 		 				.,,,,,,,,,
		1 1 1				 		 		 	1444			
		1111111111	Laboration of the same of					 						
7177		*****				 	 	 		 				
						 11111111	 	 		 				
1111	1000						 	 		 				
1111						 	 	 		 				
::::		4 5 5 5 5 5 5 5 5			1	 	 	 						

Name:	 	
Name:	 	······
Name:	 	······································
Group:	 Note:	/